



Graduated beakers, PP, raised blue scale



Highly transparent. With easily readable, raised, embossed blue scale, and stable, easy-grip handle. To preserve markings, do not clean at temperatures exceeding 60 °C. Conditionally autoclavable at 121 °C (2 bar) according to DIN EN 285. For autoclaving, we recommend the design with raised graduations (cat. nos. 440941 - 447941).

Volume ml	Divisions ml	Height mm	Ø mm	PU	Cat. No.
50	2	70	40	24	446081
100	2	80	50	24	447081
250	5	120	74	12	440081
500	10	140	92	12	441081
1000	20	181	117	6	442081
2000	50	213	152	6	443081
3000	50	242	172	6	444081
5000	100	270	204	6	445081



Graduated beakers, PP, raised scale



Highly transparent. With a raised scale and stable, easy-grip handle. Autoclavable at 121 °C (2 bar) according to DIN EN 285.

Volume ml	Divisions ml	Height mm	Ø mm	PU	Cat. No.
50	2	70	40	24	446941
100	2	80	50	24	447941
250	5	120	74	12	440941
500	10	140	92	12	441941
1000	20	181	117	6	442941
2000	50	213	152	6	443941
3000	50	242	172	6	444941
5000	100	270	204	6	445941

Measuring and transferring

Graduated beakers, stackable, PP



Highly transparent. With stable handle and easily readable, printed black scale on both sides. Therefore, the volume is equally visible for left- and right handers.

To preserve markings, cleaning at no higher than 60 °C is recommended. For autoclaving we recommend the design with raised graduations (cat. nos. 440941 - 447941).



Volume ml	Divisions ml	Height mm	Ø mm	PU	Cat. No.
250	5	115	75	12	480941
500	10	140	100	12	481941
1000	10	167	125	12	482941
2000	20	212	148	12	483941
3000	50	242	170	12	484941

Graduated beakers, stackable, coloured, PP



Graduated beakers in four different colours. Transparent, with stable handle and easily readable, printed scale on both sides. Therefore, the volume is equally visible for left- and right handers. To preserve markings, cleaning at no higher than 60 °C is recommended. For autoclaving we recommend the design with raised graduations (cat. nos. 440941 - 447941).



Volume ml	Colour	Divisions ml	Height mm	Ø mm	PU	Cat. No.
500	blue	10	140	100	12	481942
500	yellow	10	140	100	12	481943
500	red	10	140	100	12	481944
500	green	10	140	100	12	481945
500	Set: blue, yellow, red, green (1 item each)	10	140	100	1	4811111
1000	blue	10	167	125	12	482942
1000	yellow	10	167	125	12	482943
1000	red	10	167	125	12	482944
1000	green	10	167	125	12	482945
1000	Set: blue, yellow, red, green (1 item each)	10	167	125	1	4821111



Graduated beakers, SAN



Crystal clear.
With a raised scale and stable, easy-grip handle.

Volume ml	Divisions ml	Height mm	Ø mm	PU	Cat. No.
250	5	120	70	12	44091
500	10	133	91	12	44191
1000	10	170	116	6	44291
2000	20	215	150	6	44391
3000	50	242	170	6	44491



Collectors, PP or SAN



With a raised scale. Volume: 2000 ml.
With stable, easy-grip handle and white PC lid.
Diameter: 150 mm; height: 220 mm.

Description	Divisions ml	PU	Cat. No.
SAN, raised scale (Picture 1)	20	6	97891
PP, raised scale	50	6	978941
PP, raised, blue embossed scale (Picture 2)	50	6	978081
Accessories for collectors			
Lid, PC		6	97791

