

Sample preparation

Erlenmeyer flasks, PMP with PP screw cap



Transparent. Wide-mouth, can also be used with NS stoppers (not included).

Ideal for use as a receiving vessel in titrations.

Well suited for storage and cultivation of cell cultures. Far safer than glass flasks for use in incubator shakers due to the break resistance of plastic. Suitable for microwaves.

To preserve markings, cleaning at no higher than 60 °C is recommended.



Volume ml	Divisions ml	Thread GL	Neck NS	PU	Cat. No.
50	10	40	34/35	6	66695
100	20	40	34/35	6	66795
250	50	52	45/40	6	66895
500	100	52	45/40	6	66995
1000	200	52	45/40	4	67095

Erlenmeyer flasks, PP with PP screw cap



Transparent. Wide-mouth, can also be used with NS stoppers (not included).

Well suited for storage and cultivation of cell cultures. Far safer than glass flasks for use in incubator shakers due to the break resistance of plastic. Suitable for microwaves.

To preserve markings, cleaning at no higher than 60 °C is recommended.



Volume ml	Divisions ml	Thread GL	Neck NS	PU	Cat. No.
50	10	40	34/35	6	666941
100	20	40	34/35	6	667941
250	50	52	45/40	6	668941
500	100	52	45/40	6	669941
1000	200	52	45/40	4	670941

NS stoppers can be found on page 112.