



Drying rack

Back plate and trough made from PVC with drainage nozzle.
 With 75 metal pegs (length: 10 cm) with PE coating for hanging various sizes of apparatus.
 With two bore holes for simple wall mounting.
 Delivered without installation hardware.

Size mm	PU	Cat. No.
450 x 630	1	76299



Drying rack, PS

With wide draining trough and drainage nozzles.
 Rack with 72 pegs 95 x 15 mm. For drying larger apparatus, some of the pegs can be removed and the bore holes closed at the rear.
 Delivered complete with drainage tube and accessories for the wall installation.
 In addition, 11 pegs (95 x 6 mm) are included for objects having a smaller diameter, such as reagent tubes.

Description	Size mm	PU	Cat. No.
Drying rack	450 x 630	1	80213
Pegs	95 x 6	11	81213